


PATIENT MONITORING



VERMED[®]
a Graphic Controls company

Patient Monitoring

Vermed manufactures and provides dependable, high quality and easy to use consumable medical products. Our products are found throughout hospitals around the world enhancing patient care.




ECG ELECTRODES

Electrode Anatomy	4
Product Packaging	6
VersaTrode	8
BluFlex	10
Swaromed	12
Swaromed Pediatric	14
Prewired Pediatric	16
Lead-Locking	18
SilveRest Tabs	20




DEFIB PADS

Defib Pad Anatomy	22
Philips Medical	24
Physio Control	24
Zoll Medical	24
Cardiac Science	24




GROUNDING PADS (Diathermy Pad | Neutral Electrodes)

Swaroplate Non-Split	26
Swaroplate Split	26
Prewired	26
Twin-Safe	26



MEDICAL CHART PAPERS

Burdick	29	Instrumentation Labs	31	Quinton	32
Corometrics	29	Marquette GE	31	Radiometer	32
Datascope	29	Medtronic	31	Schiller	32/33
Dinamap	29	Mennen	31	Siemens	33
Edan	30	Micro Medical	31	Sonicaid (Huntleigh)	33
Electodyne	30	Mindray	31	Spacelabs	33
Fukuda Denshi	30	Miscellaneous	31	St. Johns	33
GE	30	Mortara	32	Steris	33
Geringe	30	Nihon Kohden	32	Welch Allyn	33
Hellige	30	Philips	32	Vitalograph	33
Hewlett Packard	30	Physio Control	32	Zoll	33
Humphrey	31	Protocol	32		




FACILITIES CHARTS

Helmer	35
Jewitt	35
United Electric	35
Jules Richard	35


Let Colour be Your Guide

We know that selecting the best product for every application saves you time and money - so we make it easy. Colour coding for quick reference means confidence in your selection. Every section through out the catalog has its own unique colour, allowing you to quickly and easily find exactly what you are looking for.



ULTRASOUND FILM

Sony	37
Mitsubishi	38
Vermed Ultra	39



FETAL BELTS

Abdominal Transducer Belts	41
----------------------------	----



ULTRASOUND GELS

Vermed	43
Aquasonic 100	43
Aquasonic Clear	43
SCAN Gels	44
Ultrasound Probe Covers	44
Thermasonic Gel Warmers	44
Conductive Media	44



STAINLESS SNAP

Stainless steel snap stops corrosion from occurring, eliminating artifact in the tracing. Contains a Silver / Silver Chloride sensor.

SUBSTRATE

Vermed Electrodes are one millimeter thick. Our foam electrodes provide more strength and less stretch during placement.

WET GEL

Our wet gel electrodes come in a variety of chloride levels allowing for quicker and more reliable traces.

LABEL

The dot label makes sure the snap is held in place on the electrode. It is also generally colour coordinated to help identify product lines and application.

CARBON SNAP

Our Carbon snaps are MR Conditional. It can be worn during MRI Scans, but not actively in use. Radiolucent and contains a Silver / Silver Chloride sensor.

ADHESIVE

Our adhesive material is made from a medical grade acrylic - ensuring that your electrodes stay in place. Larger electrodes allow for an increase in adhesion.

SOLID GEL

Solid gel is filled in at the depth of the foam to allow for smooth, flat application. Solid gel aids in the adhesion of the electrode and does not require clean up after use.



LINERS

Wet gel and solid gel electrodes on tape or cloth come on an opaque liner that has a thermally molded cup to protect the integrity of the gel and prevent it's dispersion on the adhesive area.

GEL OPTIONS



Solid Gel




Wet Gel

SNAP OPTIONS



Stainless Steel



Carbon Filled

MR Conditional



SUBSTRATE OPTIONS

Foam Electrodes

1mm Polyethylene Foam

Cloth Electrodes

Spunlace Non-Woven Fabric

BluFlex Electrodes


Spunlace Non-Woven Fabric




Latex Free

vermed.co.uk

RESEALABLE BULK POUCHES
25 or more electrodes.



BULK PACK (Front)
Country Specific Languages




BULK PACK (Back)
Pictogram and Country
Specific Languages.
Directions For Use (DFU's)




CASE LABELS
Product # and information
placed on case packaging

INDIVIDUAL ROLLS
Individually shrink wrapped rolls
bulk packed within the case.




INNER SLEEVES
5 individual rolls are packaged in
inner sleeves within the case.



MULTIROLL PACK
Sets of 5 or 10 rolls shrink
wrapped together.



FAN FOLD
Individually shrink wrapped
fanfold packs within a case.



ULTRASOUND FILM CASE
Individually packaged rolls of
film packed within the case.





Vermed's VersaTrode™ line of disposable ECG electrodes are the ideal choice for all adult cardiac monitoring and diagnostic applications. Each VersaTrode™ is designed to work in all procedures: Stress Test, Long Term Monitoring, Short Term Monitoring or Resting ECG.

60mm - Round - Foam Electrodes

Product #	Reference #	Per Pouch	Per Box	Per Case	Snap	Adhesive	Gel Type
6000937N	FW60	25	n/a	250	Stainless Steel	Regular	Wet Gel 
6000939N	HTFW60	25	n/a	250	Stainless Steel	High Tack	Wet Gel 

50mm - Round - Foam Electrodes

Product #	Reference #	Per Pouch	Per Box	Per Case	Snap	Adhesive	Gel Type
6000936N	FW50	50	n/a	400	Stainless Steel	Regular	Wet Gel 
6000935N	FW50R	50	n/a	400	PP Carbon Blend	Regular	Wet Gel  
6000941N	FS50	50	n/a	400	Stainless Steel	Regular	Solid Gel 
6000950N	FS50R	50	n/a	400	PP Carbon Blend	Regular	Solid Gel  

36 x 50mm - Oval - Foam Electrodes

Product #	Reference #	Per Pouch	Per Box	Per Case	Snap	Adhesive	Gel Type
6000972N	HTFW3654	25	n/a	500	4mm	High Tack	Wet Gel 
6000982N	HTFW3654	50	n/a	400	4mm	High Tack	Wet Gel  
6000940N	HTFS365	50	n/a	500	Stainless Steel	High Tack	Solid Gel 

32mm - Square - Foam Electrodes

Product #	Reference #	Per Pouch	Per Box	Per Case	Snap	Adhesive	Gel Type
6000930N	HTFW32	50	n/a	500	Stainless Steel	High Tack	Wet Gel 
6001023N	EFS32	50	n/a	500	Stainless Steel	High Tack	Solid Gel 
6001126N	FS32/R	50	n/a	500	PP Carbon Blend	High Tack	Solid Gel  
6001200N	FW32323RBP	50	n/a	400	Stainless Steel	High Tack	Solid Gel  

32mm - Teardrop - Foam Electrodes

Product #	Reference #	Per Pouch	Per Box	Per Case	Snap	Adhesive	Gel Type
6000931N	FW32P	50	n/a	500	Stainless Steel	Regular	Wet Gel 



FW60 | HTFW60
60mm



FW50R | FS50R
50mm Radiolucent



HTFS365 | HTFW3654
36 x 50mm



HTFW32 | ESF32 | FS32R | FW32
32mm (Radiolucent Option)

FW32P
32mm

FEATURES

- Adult foam wet and solid gel electrodes
- Hypoallergenic adhesive
- High Tack option offers stronger adhesive
- Good for use throughout the hospital





BluFlex electrodes are made of a very soft and flexible material, allowing for maximum comfort for your patient. Made of Hydrophobic Spun Lace creating a soft, flexible material optimizing patient comfort during extended wear applications. Our BluFlex material is fluid resistant and breathable, making it comparable to tape and cloth.

60mm - Round - Nonwoven Cloth Electrodes

Product #	Reference #	Per Pouch	Per Box	Per Case	Snap	Adhesive	Gel Type	
6000036N	1033	25	250	3000	Stainless Steel	Regular	Wet Gel	💧
6000636N	1033	30	300	3600	Stainless Steel	Regular	Wet Gel	💧
6000929N	BW60	25	250	3000	Stainless Steel	Regular	Wet Gel	💧

55mm - Round - Nonwoven Cloth Electrodes

Product #	Reference #	Per Pouch	Per Box	Per Case	Snap	Adhesive	Gel Type	
6000032N	1032	25	250	3000	Stainless Steel	Regular	Wet Gel	💧
6000432N	1032	50	400	4800	Stainless Steel	Regular	Wet Gel	💧
6000658N	1032	30	300	3600	Stainless Steel	Regular	Wet Gel	💧

50mm - Round - Nonwoven Cloth Electrodes

Product #	Reference #	Per Pouch	Per Box	Per Case	Snap	Adhesive	Gel Type	
6000051N	1039	25	250	3000	Stainless Steel	Regular	Wet Gel	💧
6000430N	1039	50	400	4800	Stainless Steel	Regular	Wet Gel	💧
6000932N	BW50	50	400	4800	Stainless Steel	Regular	Wet Gel	💧

40mm - Round - Nonwoven Cloth Electrodes

Product #	Reference #	Per Pouch	Per Box	Per Case	Snap	Adhesive	Gel Type	
6000016N	1037D	25	250	3000	Stainless Steel	Regular	Wet Gel	💧
6001021N	1037D	50	400	4800	Stainless Steel	Regular	Wet Gel	💧

32mm - Teardrop - Nonwoven Cloth Electrodes

Product #	Reference #	Per Pouch	Per Box	Per Case	Snap	Adhesive	Gel Type	
6000041N	1035	25	250	3000	Stainless Steel	Regular	Wet Gel	💧
6000934N	BW32P	50	500	6000	Stainless Steel	Regular	Wet Gel	💧



FEATURES

- Adult spun lace wet gel electrode for long term monitoring
- Tape, foam and cloth compatible
- Hypoallergenic adhesive
- Wet Gel offers higher conductivity and lower impedance
- 4 to 5 days of extended wear
- Good for use in most longer term applications






FEATURES

Swaromed ECG Electrodes are pre-gelled (single-use) and have a Silver/Silver Chloride sensor. All adhesive foams and types are equipped with a hypoallergenic medical-grade adhesive that is non-irritating to the skin.

- Hypoallergenic adhesive
- Silver/Silver Chloride sensor
- Quick recovery gel
- Strict quality control during production
- Defibrillation resistant



Swaromed™ Adult General Purpose Electrodes

60mm - Round - Foam Electrodes

Product #	Reference #	Per Pouch	Per Box	Per Case	Snap	Adhesive	Gel Type
6000937N	1001	25	250	3000	Stainless Steel	Regular	Wet Gel
6001122N	4001	25	250	3000	Stainless Steel	Regular	Solid Gel
6000004N	1002	25	250	3000	Stainless Steel	High Tack	Wet Gel
6000003N	4002	25	250	3000	Stainless Steel	High Tack	Solid Gel

50mm - Round - Foam Electrodes

Product #	Reference #	Per Pouch	Per Box	Per Case	Snap	Adhesive	Gel Type
6000017N	1009	25	250	3000	Stainless Steel	Regular	Wet Gel
6000018N	1009	50	400	4800	Stainless Steel	Regular	Wet Gel
6000610N	1009	30	300	3600	Stainless Steel	Regular	Wet Gel
6000388N	1009A	25	250	3000	Stainless Steel	Regular	Wet Gel
6000582N	4009	50	400	4800	Stainless Steel	Regular	Solid Gel
6000020N	1010	25	250	3000	PP Carbon Blend	Regular	Wet Gel
6000436N	1010	50	400	4800	PP Carbon Blend	Regular	Wet Gel
6000611N	1010	30	300	3600	PP Carbon Blend	Regular	Wet Gel
6000981N	4010	50	400	4800	PP Carbon Blend	Regular	Solid Gel
6000665N	1012	50	400	4800	PP Carbon Blend	High Tack	Wet Gel
6000025N	1019	25	250	3000	Stainless Steel	High Tack	Wet Gel
6000341N	1019	50	400	4800	Stainless Steel	High Tack	Wet Gel
6000605N	1019	30	300	3600	Stainless Steel	High Tack	Wet Gel
6000616N	4019	50	400	4800	Stainless Steel	High Tack	Solid Gel

55mm - Round - Foam Electrodes

Product #	Reference #	Per Pouch	Per Box	Per Case	Snap	Adhesive	Gel Type
6000028N	1023	25	250	3000	Stainless Steel	Regular	Wet Gel
6000059N	1023	50	400	4800	Stainless Steel	Regular	Wet Gel
6000607N	1023	30	300	3600	Stainless Steel	Regular	Wet Gel
6000397N	1022	50	400	4800	Stainless Steel	High Tack	Wet Gel
6000652N	4022	50	400	4800	Stainless Steel	High Tack	Solid Gel

36 x 50mm - Oval - Foam Electrodes

Product #	Reference #	Per Pouch	Per Box	Per Case	Snap	Adhesive	Gel Type
6000045N	1036	50	500	6000	Stainless Steel	High Tack	Wet Gel
6000608N	1036	30	300	3600	Stainless Steel	High Tack	Wet Gel
6000679N	1036A	50	500	6000	Stainless Steel	High Tack	Wet Gel
6000578N	4036	50	500	6000	Stainless Steel	High Tack	Solid Gel
6000398N	1040	50	500	6000	PP Carbon Blend	High Tack	Wet Gel
6000069N	1084	25	250	3000	Stainless Steel	High Tack	Wet Gel
6000609N	1084	30	300	3600	Stainless Steel	High Tack	Wet Gel
6000438N	1084FT	30	300	3600	Stainless Steel	High Tack	Wet Gel





40mm - Round - Foam Electrodes

Product #	Reference #	Per Pouch	Per Box	Per Case	Snap	Adhesive	Gel Type
6000446N	1007D	50	400	4800	PP Carbon Blend	Regular	Wet Gel 
6000495N	1008D	50	400	4800	Stainless Steel	Regular	Wet Gel 
6000458N	4008D	50	500	6000	Stainless Steel	Regular	Solid Gel 
6000614N	4008D	30	300	3600	Stainless Steel	Regular	Solid Gel 
6000973N	4008D	50	400	4800	Stainless Steel	Regular	Solid Gel 

32mm - Teardrop - Foam Electrodes

Product #	Reference #	Per Pouch	Per Box	Per Case	Snap	Adhesive	Gel Type
6000063N	1057	25	250	300	Stainless Steel	Regular	Wet Gel 
6000612N	1057	30	300	3600	Stainless Steel	Regular	Wet Gel 

32mm - Square - Foam Electrodes

Product #	Reference #	Per Pouch	Per Box	Per Case	Snap	Adhesive	Gel Type
6000044N	1060	50	500	6000	PP Carbon Blend	High Tack	Wet Gel  
6000588N	4060	50	500	6000	PP Carbon Blend	High Tack	Solid Gel  
6000587N	1088	50	400	4800	PP Carbon Blend	High Tack	Wet Gel  
6000065N	1066	50	500	6000	Stainless Steel	High Tack	Wet Gel 
6000534N	1066A	50	500	6000	Stainless Steel	High Tack	Wet Gel  
6000581N	4066	50	500	6000	Stainless Steel	High Tack	Solid Gel 
6000068N	1083	25	250	300	Stainless Steel	High Tack	Wet Gel  
6000570N	1083	50	400	4800	Stainless Steel	High Tack	Wet Gel  



1007D | 1008D | 4008D - Diddle
40mm Pediatric (Radiolucent Option)



1057
32mm Pediatric



1066 | 4066
32mm Pediatric



1060 | 4060
32mm Radiolucent



1088 | 1083
32mm Banana Plug

FEATURES

- Smaller electrodes for easier placement
- Available in wet or solid gel options
- Doodle on label helps smaller patients feel more comfortable
- Available with stainless steel or radiotranslucent PP Carbon snaps
- Banana plug adapters available

Abrader Cap



Banana Plug



Latex-free



Wet Gel



Solid Gel



MR Conditional






Vermed offers a complete line of pre-wired electrodes for use in the PICU/NICU and other neonatal applications. Available in various sizes as well as metallic or radiolucent lead wires.

Pre-Wired Neonatal ECG

Product #	Description	Material	Per Bag	Per Box		
4009446C	23 x 23 x 1mm Square	Cloth	3	150	Solid Gel	MR Conditional
4009448C	23 x 23 x 1mm Square 1.5mm connector	Cloth	3	150	Solid Gel	
4009449C	23 x 23 x 1mm Square 4mm connector	Cloth	3	150	Solid Gel	
4009450C	23 x 23 x 1mm Oval	Cloth	3	150	Solid Gel	MR Conditional



FEATURES

Pre-Wired Electrodes:

We offer a complete line of pre-wired electrodes for use in PICU/NICU and other neonatal applications. Available in various sizes as well as metallic or radiolucent lead wires.

MR Conditional




Latex-free



Solid Gel





Don't let wires get in the way. LeadLocking™ electrodes save time and money by not having to tape down wires or replace displaced electrodes. Ensure patient comfort and accuracy by preventing motion artifact.

6001246N - 55mm x 40mm - Rectangular foam
LeadLocking™ Electrode with metal snap

Product #	Per Strip	Per Pouch	Per Box	Per Case	Solid Gel
6001246N	n/a	30	n/a	300	

6001252N - 55mm x 40mm - Rectangular foam
LeadLocking™ Electrode with metal snap

Product #	Per Strip	Per Pouch	Per Box	Per Case	Wet Gel
6001252N	n/a	30	n/a	300	

FEATURES

- Wire restraint electrode safely secures leadwires preventing artifact associated with movement
- Specially designed electrode shape helps ease of correct placement
- Stainless steel snap minimizes the amount of noise artifact due to lower incidence of corrosion compared to nickel plated brass used by competitors



6001246N
Foam | Solid Gel



6001252N
Foam | Wet Gel



SilveRest™ Tab Electrodes

Diagnostic Resting EKG Tab Electrode

Vermed SilveRest™ Tab Electrodes for resting EKG provide excellent adhesion and quality traces. All of our tab electrodes meet or exceed the AAMI standard for disposable pre-gelled EKG electrodes.

4007847C - 26.4mm x 22.5mm - Ag/AgCL

- Safe, smooth and comfortable adhesion that conforms to body shape.
- Easy application and use helps provides a more reliable trace .

Product #	Per Strip	Per Pouch	Per Box	Per Case	Solid Gel
4004847C	10	100	2500	5000	

A10201 - 1.34" x .90" / 34.036mm x 22.86mm - C/AgCl

- Proprietary low residue gel for easy clean-up.
- Larger clip attachment area uses small or large clips without compromise.

Product #	Per Strip	Per Pouch	Per Box	Per Case	Solid Gel
A10201	10	100	1000	4000	

A10023 - 1.125" x 1.438" / 28.575mm x 36.525mm - Ag/AgCL

- Proprietary adhesive solid gel for optimal adhesion.
- Larger surface area contact which reduces impedance and provides superior adhesion.
- Repositionable and packaged in resealable pouches

Product #	Per Strip	Per Pouch	Per Box	Per Case	Solid Gel
A10023	10	100	1000	4000	

FEATURES

- Conductive adhesive gel
- Patented traction liner prevents tabs from sliding off the liner - saving time and money
- Unique liner card ensures packs do not brick together for ease of removal and application
- Pliable substrate conforms to the body
- Adhesive solid gel requires minimal clean-up
- Latex-free

Latex-free



Solid Gel




4007847

Tab | Solid Gel
Silver /Silver



A10201


Tab | Solid Gel
Silver / Carbon



A10023

Regal Tab | Solid Gel
Silver / Silver





60" LEAD WIRES

Longest in the market. Provides more flexibility when placing the defibrillator.

Colour CODED PADS

Anterior electrode is blue foam and the posterior electrode is white foam. This helps train staff to look for colours to aid in pad placement.



PAD SIZE

Our pads are smaller than competitor pads, allowing for easier placement on all body shapes and sizes.

CONDUCTOR PLATE

Our aggressive adhesive conductive solid hydrogel reduces the need for a large foam boarder allowing for a larger conductor plate than the competition.



INSTRUCTIONS FOR USE

User instructions and a check list is on the coloured release liner. Coloured liner does not become a slipping hazard for the staff on the treatment floor.

Instructions for use:

1. Remove electrodes from peel pouch. ☒
2. Ensure skin site is clean and dry. ☐
3. Remove electrodes from release liner. ☐
4. Apply electrodes firmly to patient's skin. ☐
5. Connect electrodes to defibrillator or patient cable. ☐
6. Connect ECG leads for demand pacing. ☐

PEEL
HERE
ONLY

PEEL
HERE
ONLY

SHELF LIFE

Vermed Defib Pads have a 36 month shelf life (3 years). That's approximately 50% longer than our leading competitors.

POUCH OPENING

Vermed Defib Pads have a traditional notched tear opening, or an easy to peel alternative, providing ease of use for clinicians.

Colour CODED PACKAGING

Eliminate mistakes. Our labels are consistent with Medtronic in red, Philips in green, Zoll in yellow, HP in orange and Anderson is blue. Our Adult Pouch background is white and our Pediatric pouch background is light blue.

CONNECTOR STANDARDIZATION

Ask about our Connector Standardization Program to help standardize all your connectors across an entire hospital. This allows for one defib pad to be used in all departments, and reduces errors in choosing the proper pad.



Defibrillation Electrodes

Defib Pads for All Departments of the Hospital

Vermed's Defib Pads were designed with your needs in mind.
Our Defib Pad line is not just different - it is a step above the other brands

PHILIPS MEDICAL

Product #	For Defibrillator Models	Pad & Leadwires
D311	XLT, MRX, HeartStart XL, MRX ALS, HP Codemaster XL, HeartStart FR2+, HeartStart Forerunner	Radiotransparent - Invisible
D309	XLT, MRX, HeartStart XL, MRX ALS, HP Codemaster XL, HeartStart FR2+, HeartStart Forerunner	Radiolucent - Slight Shadow
D317	XLT, MRX, HeartStart XL, MRX ALS, HP Codemaster XL	Radiolucent - Slight Shadow Leads Out Package
D321	XLT, MRX, HeartStart XL, MRX ALS, HP Codemaster XL	Pediatric Radiolucent - Slight Shadow

PHYSIO CONTROL

Product #	For Defibrillator Models	Pad & Leadwires
D304	Lifepak 9P, 15, 10C, 20, 12, 500, 500 DPS, 1000	Radiotransparent - Invisible
D313	Lifepak 9P, 15, 10C, 20, 12, 500, 500 DPS, 1000	Radiolucent - Slight Shadow
D315	Lifepak 9P, 15, 10C, 20, 12, 500, 500 DPS, 1000	Radiolucent - Slight Shadow Leads Out Package
D319	Lifepak 9P, 15, 10C, 20, 12	Pediatric Radiolucent - Slight Shadow

ZOLL MEDICAL

Product #	For Defibrillator Models	Pad & Leadwires
D303	M-Series, E-Series, R-Series, X-Series	Radiotransparent - Invisible
D312	M-Series, E-Series, R-Series, X-Series	Radiolucent - Slight Shadow
D314	M-Series, E-Series, R-Series, X-Series	Radiolucent - Slight Shadow Leads Out Package
D318	M-Series, E-Series, R-Series, X-Series	Pediatric Radiolucent - Slight Shadow

CARDIAC SCIENCE

Product #	For Defibrillator Models	
D326	PowerHeart AED G3 Plus, PowerHeart AED G3 Pro	Can be used in place of Cardiac Science Part #9131-001


FEATURES

- 36 Month Shelf Life
- Proprietary Adhesive Gel
- Larger Conductor Plate
- 60" Lead Wires
- Smaller Pad Size
- Connector Standardization


Latex-free Solid Gel




PHILIPS MEDICAL Connector Style




PHYSIO CONTROL Connector Style



ZOLL MEDICAL Connector Style



CARDIAC SCIENCE Connector Style



The OEM trademarks identified herein are the trademarks of the respective OEMs, and not Vermed. Vermed disclaims any affiliation, connection or association between its products and those of the respective OEMs.



Swaroplate Non-Split

Cross Ref.	Use	Size (Length/Width)	Adhesive Surface	Conductive Surface	Impedance	Temp. Rise
1201H	Adults	215 x 104mm	207 cm ²	134 cm ²	1 Ohm	< 3°C
1202H	Children	144 x 90mm	115 cm ²	72 cm ²	2 Ohm	< 5°C/500mA
1203H	Infants & Prematures	165 x 74mm	66 cm ²	35 cm ²	4 Ohm	< 5°C/500mA
1210H	Adults & Children	156 x 126mm	144 cm ²	105 cm ²	1.5 Ohm	< 4°C

Swaroplate Split

Cross Ref.	Use	Size (Length/Width)	Adhesive Surface	Conductive Surface	Impedance	Temp. Rise
1211H	Adults	201 x 106mm	200 cm ²	144 cm ²	1 Ohm	< 3°C
1212H	Adults	206 x 122mm	220 cm ²	160 cm ²	1 Ohm	< 3°C
1213H	Adults	230 x 123mm	232 cm ²	169 cm ²	1 Ohm	< 3°C
1222H	Adults & Children	154 x 88mm	124 cm ²	71 cm ²	2 Ohm	< 5°C
1226H	Adults & Children	164 x 124mm	105 cm ²	71 cm ²	2 Ohm	< 5°C
1223H	Children	165 x 74mm	66 cm ²	32 cm ²	4 Ohm	< 5°C
1220H	Infants & Prematures	160 x 122mm	148 cm ²	105 cm ²	1.5 Ohm	< 4°C

Pre-Wired

Cross Ref.	Use	Size (Length/Width)	Adhesive Surface	Conductive Surface	Base Pad	Connector
1251H	Adults	201 x 106mm	200 cm ²	144 cm ²	1211	Valleylab REM
1255H	Adults	210 x 127mm	220 cm ²	160 cm ²	1212	Valleylab REM
1261H	Adults	215 x 104mm	232 cm ²	169 cm ²	1201	Valleylab NON-REM
1252H	Children	154 x 88mm	124 cm ²	71 cm ²	1222	Valleylab REM
1257H	Children	164 x 124mm	105 cm ²	71 cm ²	1226	Valleylab REM
1253H	Infants & Prematures	165 x 74mm	66 cm ²	32 cm ²	1223	Valleylab REM

Twin-Safe

Cross Ref.	Use	Size (Length/Width)	Adhesive Surface	Conductive Surface	Cable Length	Connector
1280H	Adults	176 x 122mm	169 cm ²	110 cm ²	n/a	n/a
1285H	Adults	176 x 122mm	169 cm ²	110 cm ²	3m	Valleylab REM
1286H	Adults	176 x 122mm	169 cm ²	110 cm ²	3m	Jack Connector



SWAROPLATE SPLIT
Grounding Pads



TWIN-SAFE
Grounding Pads

GROUNDING PADS

From our unique TwinSafe design to our standard Swaroplate models, Vermed Grounding (Diathermy) Pads are compatible with many RF generators. Pads are produced with the highest quality materials and designed to ensure complete contact during HF surgery.



Burdick

Product #	Cross Ref.	Description	Channels	UOM	Per Case
32017246	BDK 007868	8.438" x 183'	Single Channel	Pad	10 Pads
32017197	BDK 007966	8.438" x 183' 200sh/pd	Single Channel	Pad	10 Pads
2018511C	BDK 007958	50mm x 150'	Single Channel	Roll	10 Rolls
31334595	BDK 007983	8.465" x 183' 200sh/pd	Single Channel	Pad	10 Pads
32916653	BDK 716-0237-00	8.465" x 183' 200sh/pd	Single Channel	Pad	10 Pads
32017188	BDK 716-0234-00	8.5" x 183' 200sh/pd	Multi Channel	Pad	10 Pads
32017526	BDK 007988 (716-0239-00)	215mm x 137.5' 300sh/pd	Multi Channel	Pad	10 Pads
32017526	BDK 716-0239-00	215mm x 137.5' 300sh/pd	Multi Channel	Pad	10 Pads
2026689C	BDK 007985	210mm x 60m 200sh/pd	Multi Channel	Pad	10 Pads

Corometrics

Product #	Cross Ref.	Description	Channels	UOM	Per Case
30597226	CMS 4305	152mm x 47' with Labels	Fetal	Pad	40 Pads
30767589	CMS 4305 BAO	152mm x 47' 160pg	Fetal	Pad	40 Pads
2003454C	CMS 4483	140mm X 130'	Fetal	Pad	12 Pads
4001556C	CMS 4305DAO 25 yr	152mm x 47'	Fetal	Pad	40 Pads

Datascope


Product #	Cross Ref.	Description	Channels	UOM	Per Case
30015682	0683-00-0422-02 (SPA 307438)	50mm x 87'	Single Channel	Roll	50 Rolls

Dinamap

Product #	Cross Ref.	Description	Channels	UOM	Per Case
7G32018145	DM 89100	2.25" x 80' No Print	n/a	Roll	10 Rolls

MEDICAL CHART PAPER

VERMED is the ORIGINAL Chart Paper manufacturer, offering the widest assortment of chart papers for all departments in the hospital - from cardiology to the blood bank. We can consolidate the greatest number of suppliers for you – saving time and money. Vermed manufactures in the UK and USA to strict OEM specifications. We offer the exact same grade of paper, dimensions and archival qualities as the Original Equipment Manufacturers. In fact, we make charts for many OEMs, Offering the best VALUE product on the market - OEM quality at a significant savings.



Edan

Product #	Cross Ref.	Description	Channels	UOM	Per Case
4007378C	EDN SE-12/1200	210mm x 295mm	Multi Channel	Pad	1 Pad

Electrodyne

Product #	Cross Ref.	Description	Channels	UOM	Per Case
32017187	ELD K140-1000	50mm x 89' Black Grid	Single Channel	Roll	100 Rolls
10669183	ELD K140-1010	50mm x 89' Red Grid	Single Channel	Roll	50 Rolls
10784107	ELD K140-1020	50mm x 89' Green Grid	Single Channel	Roll	50 Rolls

Fukuda Denshi

Product #	Cross Ref.	Description	Channels	UOM	Per Case
31052726	FKD 050F30-H02	50mm x 98'	Single Channel	Pad	120 Pads
2016931C	FKD OP-119TE	63mm x 97'	Single Channel	Roll	20 Rolls
4008207C	FKD OP-219TE	63mm x 99' 302sh/pd	Multi Channel	Pad	12 Pads
2025248C	FKD OP-124TE	49.5mm x 67' 272sh/pd	Multi Channel	Pad	10 Pads
2001532C	FKD OP-69TE	209mm x 88'	Multi Channel	Roll	10 Rolls

GE

Product #	Cross Ref.	Description	Channels	UOM	Per Case
4000936C	GE 2030887-001	80mm x 82' 280sh/pd	Single Channel	Pad	10 Pads
32018479	GE 74205-HEL	49.5mm x 84' No Print	n/a	Roll	20 Rolls
4001144C	GE 2036970-001 MAC 800	110mm x28m 200sh/pad	Multi Channel	Pad	10 Pads
4000936C	GE 2030887-001 MAC 400/600	79.5mm x 25.2m	Multi Channel	Pad	10 Pads

Getinge

Product #	Cross Ref.	Description	Channels	UOM	Per Case
32018706	GTG 61301600005	2.25" x 107"(can order 8)	n/a	Roll	80 Rolls
31449534	MRC GTG 61301601272	2.36 x 70' No Print	n/a	Roll	10 Rolls

Hellige

Product #	Cross Ref.	Description	Channels	UOM	Per Case
2022033C	HEB 22616605	210mm x 45m 150sh/pd	Multi Channel	Pad	10 Pads
4000468C	HEB 22616611 25 yr	210mm x 45m 25yr	Multi Channel	Pad	1 Pad

Hewlett Packard

Product #	Cross Ref.	Description	Channels	UOM	Per Case
10005156	HP 9270-0484	151.6mm x 50'	Fetal	Pad	40 Pads
10632900	HP 40457A (can order 10)	50mm x 100'	Single Channel	Roll	50 Rolls
2000188C	HP 9270-0485	151.6mm x 50'	Fetal	Pad	40 Pads
20060622	HP 9270-0980 (can order 1)	50mm x 100'	Single Channel	Roll	50 Rolls
2000036C	HP 40477A	50mm x 89' No Print	Single Channel	Roll	10 Rolls
30748696	HP M1910A	5.894" x 49'	Fetal	Pad	40 Pads
30768678	HP M1707A	8.49" x 183' 200sh/pd	Multi Channel	Pad	8 Pads
4006120C	HP M1709A	209.5mm x 197'	Multi Channel	Pad	8 Pads
31104162	HP M2481A	8.484" x 183'	Single Channel	Pad	8 Pads
32016831	HP M1911A	150mm x 49'	Fetal	Pad	40 Pads
32021183	HP M1911A 25yr	150mm x 49' 25yr	Fetal	Pad	40 Pads

Hewlett Packard (Cont.)

Product #	Cross Ref.	Description	Channels	UOM	Per Case
32018398	HP M1914A	150mm x 49'	Fetal	Pad	60 Pads
4000700	HP 3708A	210mm x 30m 100sh/pd	Multi Channel	Pad	1 Pad

Humphrey

Product #	Cross Ref.	Description	Channels	UOM	Per Case
2010895C	HUM 24433	216mm x 30m	n/a	Roll	1 Roll

Instrumentation Labs

Product #	Cross Ref.	Description	Channels	UOM	Per Case
32018666	ILS GEM 5508	59.5mm x 90' No Print	n/a	Roll	5 Rolls

Marquette | GE

Product #	Cross Ref.	Description	Channels	UOM	Per Case
113822	MQE 9402-018	4.232" x 150'	Single Channel	Roll	32 Rolls
2000060C	MQE 9402-020	8.441" x 275' 300sh/pd	Multi Channel	Pad	1 Pad
2000151C	MQE 9402-021	8.234" x 275' 300sh/pd	Multi Channel	Pad	1 Pad
4006477C	MQE 9402-023	107.5mm x 75'	Multi Channel	Roll	2 Rolls
4004771C	MQE 9402-024	8.441" x 275' 300sh/pd	Multi Channel	Pad	8 Pads
31002226	MQE 9402-024 GRN	8.443" x 275' 300sh/pd	Multi Channel	Pad	8 Pads
31002200	MQE 9402-029	8.443" x 275'	Single Channel	Pad	8 Pads
31019568	MQE 9402-30	8.443" x 275' 300sh/pd	Multi Channel	Pad	8 Pads
32017205	MQE 9402-034 ARCHIVAL	8.433" x 275'	Multi Channel	Pad	12 Pads
2022282C	MQE 9402-037 ARCHIVAL	8.433" x 275'	Multi Channel	Pad	8 Pads
30781747	MQE 9402-046 (can order 10)	50mm x 89' No Print	Single Channel	Roll	50 Rolls
32017196	MQE 9402-053	8.5" x 11' 500sh/pd	Multi Channel	Pad	5 Pads
32016997	MQE 9402-061	8.433" x 139' 152sh/pd	Multi Channel	Pad	16 Pads

Medtronic

Product #	Cross Ref.	Description	Channels	UOM	Per Case
32020245	MEN 6092	4.33" x 100' 204ah/pd	Multi Channel	Pad	6 Pads

Mennen

Product #	Cross Ref.	Description	Channels	UOM	Per Case
30707080	MEG800-040-160	2.475" x 200'	Single Channel	Roll	50 Rolls

Micro Medical

Product #	Cross Ref.	Description	Channels	UOM	Per Case
0103931C	MI PSA1600	112mm x 8m	n/a	Roll	5 Rolls

Mindray

Product #	Cross Ref.	Description	Channels	UOM	Per Case
4001129C	MRY 0683-00-0422-02 TR Paper	50mm x 18.75m	Single Channel	Roll	5 Rolls

Miscellaneous

Product #	Cross Ref.	Description	Channels	UOM	Per Case
32016732	MED 705225 (LifeSigns 5200) (Can order 10)	57mm x 28'	Single Channel	Roll	50 Rolls

Mortara

Product #	Cross Ref.	Description	Channels	UOM	Per Case
30758273	MRN 9100-007	8.465" x 458' 500sh/pd	Multi Channel	Pad	1 Pad
32020410	MRN 9100-025-50	215mm x 92' 100sh/pd	Multi Channel	Pad	10 Pads
4000482C	MRN 9100-026-50	8.433" x 62.72m 246/sh/pd	Multi Channel	Pad	1 Pad
0104415C	MRN 9100-026-51	8.250" x 75m 250sh/pd	Multi Channel	Pad	5 Pads
4003635C	MRN 9100-028-50	108mm x 140mm 200sh/pd	Multi Channel	Pad	15 Pads

Nihon Kohden

Product #	Cross Ref.	Description	Channels	UOM	Per Case
2000064C	NKN RQS50-3	50mm x 30m	Single Channel	Roll	20 Rolls
2010078C	NKN RQS63-3	63mm x 30m	Single Channel	Roll	5 Rolls
2026452C	NKN FQW110-2-140	110mm x 21m	Multi Channel	Pad	1 Pad
2026712C	NKN FQW210-30-140	210mm x 30m	Single Channel	Pad	10 Pads

Philips

Product #	Cross Ref.	Description	Channels	UOM	Per Case
32018225	PHM 989803137011	8.464" x 11" 250sh/pd	Multi Channel	Pad	12 Pads
4000650C	PHM 989803138171	75mm x 25m	Single Channel	Roll	5 Rolls
4000651C	PHM 989803102321	215mm x 84m 300sh/pd	Multi Channel	Pad	5 Pads

Physio Control

Product #	Cross Ref.	Description	Channels	UOM	Per Case
2023274C	PHY 805319-03	106.5mm x 75'	Single Channel	Roll	2 Rolls
30718970	PHY 804700 (can order 5)	50mm x 100'	Single Channel	Roll	50 Rolls

Protocol

Product #	Cross Ref.	Description	Channels	UOM	Per Case
32017249	PRT 008-0040-98 (Can order 40)	60mm x 87' No Print	Single Channel	Roll	100 Rolls

Quinton

Product #	Cross Ref.	Description	Channels	UOM	Per Case
0105967C	QUT 036868-001	210mm x 300mm 200sh/pd	Multi Channel	Pad	10 Pads
30561354	QUT 16211-002	8.5" x 183' 200sh/pd	Multi Channel	Pad	20 Pads
30729225	QUT 18416-001	8.5" x 183' 200sh/pd	Multi Channel	Pad	20 Pads
01083203	QUT 15022-001	8.484" x 200' 200sh/pd	Multi Channel	Pad	10 Pads
2011421C	QUT 30146-001	210mm x 60m 200sh/pd	Multi Channel	Pad	1 Pad

Radiometer

Product #	Cross Ref.	Description	Channels	UOM	Per Case
2025244C	RMR 984-070	111.5mm x 142'	Single Channel	Roll	5 Rolls

Schiller

Product #	Cross Ref.	Description	Channels	UOM	Per Case
2000319C	SCH 2.157005	145mm x 115'	Single Channel	Pad	25 Pads
30999067	SCH 2.157007	70mm x 59' 200sh/pd	Single Channel	Pad	36 Pads
2000320C	SCH 2.157009	210mm x 100.8m	Multi Channel	Pad	12 Pads
2022430C	SCH 2.157014	90mm x 36m	Single Channel	Pad	1 Pad
2023563C	SCH 2.157016	90mm x 23m	Single Channel	Pad	12 Pads

Schiller (Cont.)

Product #	Cross Ref.	Description	Channels	UOM	Per Case
2023699C	SCH 2.157017	8.243" x 56m	Multi Channel	Pad	1 Pad
32016808	SCH 2.157023	210mm x 69' 150sh/pd	Multi Channel	Pad	20 Pads
0101533C	SCH 2.157025	8.243" x 48m	Multi Channel	Pad	10 Pads
0103082C	SCH 2.157026	80mm x 21m	Multi Channel	Pad	1 Pad

Siemens

Product #	Cross Ref.	Description	Channels	UOM	Per Case
20371540	SIS 43-55-060	60mm x 150'	Single Channel	Roll	50 Rolls
32018703	SIS 4398	4.37" x 80'	n/a	Roll	5 Rolls
4009491C	SIS 4698	4.37" x 80' No Print	n/a	Roll	5 Rolls
32018358	SIS 5773 (Clinitek 50/100/500)	2.24" x 79' No Print	n/a	Roll	5 Rolls
32019289	SIS 101404	2.250" x 150' No Print	n/a	Roll	10 Rolls

Sonicaid (Huntleigh)

Product #	Cross Ref.	Description	Channels	UOM	Per Case
2000089C	SON 8400-8003	142.5mm x 45m	Fetal	Pad	1 Pad
4003692C	SON 8400-8003 25yr	142.5mm x 45m 25yr	Fetal	Pad	1 Pad

Spacelabs

Product #	Cross Ref.	Description	Channels	UOM	Per Case
30015682	SPA 307438 (Can order 10)	50mm x 90'	Single Channel	Roll	50 Rolls
30540499	SPA 307568	120mm x 150' 305sh/pd	Multi Channel	Pad	20 Pads
30763711	SPA 006-0196-00	50mm x 90'	Single Channel	Pad	100 Pads
4004322C	SPA 006-0210-02	57mm x 10m	n/a	Roll	12 Rolls

St. Johns

Product #	Cross Ref.	Description	Channels	UOM	Per Case
4000461C	SJM 60003605-001	210mm x 28m	n/a	Pad	12 Pads

Steris

Product #	Cross Ref.	Description	Channels	UOM	Per Case
32016258	MRC AMS P129359-008	2.28" x 75' No Print (can order 5)	n/a	Roll	100 Rolls
32008832	AMS P764330-673	American Sterilizer	n/a	Roll	24 Rolls

Welch Allyn

Product #	Cross Ref.	Description	Channels	UOM	Per Case
4000944C	WA 94016-0000	8.243" x 56m 200sh/pd	Single Channel	Pad	5 Pads
4000505C	WA 94018-0000	8.243" x 56m 200sh/pd	Single Channel	Pad	1 Pads

Vitalograph

Product #	Cross Ref.	Description	Channels	UOM	Per Case
4007874C	VTO 44058	111mm x 35m	n/a	Roll	5 Rolls
30868940	VTO 20200	7.8125 x 8.594"	Single Channel	Pad	15 Pads

Zoll

Product #	Cross Ref.	Description	Channels	UOM	Per Case
2026560C	ZL 8000-0300	90mm x 18m 200sh/pd	Single Channel	Pad	1 Pad
4009107C	ZL 8000-000901 (X-Series)	79mm x 27.5m	Multi Channel	Roll	20 Rolls



Helmer Blood Bank Refrigeration Charts

Product #	Cross Ref.	Description	Charts / Box	Minimum Order
0105872C	HER 220419 7 Day	4.5"	52	3 Boxes
0098882C	HER 220366 7 Day	4.5"	52	3 Boxes
0101104C	HER 220064 7 Day	4.5"	52	3 Boxes
0091359C	HER 220273 7 Day	4.5"	52	1 Box

Jewitt

Product #	Cross Ref.	Description	Charts / Box	Minimum Order
0099041C	JT RDR021 - Circle Chart	7.438"	52	3 Boxes
0097630C	JT RDR020 - Circle Chart	7.438"	52	3 Boxes
0098007C	JT RDR019 - Circle Chart	7.438"	52	3 Boxes
0097831C	JT RDR068 - Circle Chart	7.438"	52	3 Boxes
0097830C	JT RDR069 - Circle Chart	7.438"	52	3 Boxes

United Electric

Product #	Cross Ref.	Description	Charts / Box	Minimum Order
0084845C	UD 6282-169 - Circle Chart	6"	50	1 Box
0087720C	UD 6282-224 - Circle Chart	6"	50	3 Boxes
0096304C	UD 6282-176 - Circle Chart	6"	50	3 Boxes
0087703C	UD 6282-240 - Circle Chart	6"	50	3 Boxes
0087561C	UD 6282-161 - Circle Chart	6"	50	3 Boxes
0096344C	UD 6282-301 - Circle Chart	6"	50	3 Boxes


Jules Richard

Product #	Cross Ref.	Description	Charts / Box	Minimum Order
0075593C	RIC D32170 - Circle Chart	4.937"	100	1 Box
0084442C	RIC D32870 - Circle Chart	4.937"	100	1 Box
0100906C	RIC D33070 - Circle Chart	4.937"	100	1 Box
0078013C	RIC D31740 - Circle Chart	4.937"	100	1 Box
0089870C	RIC D32840 - Circle Chart	4.937"	100	3 Boxes
0101142C	RIC D32340 - Circle Chart	4.937"	100	3 Boxes


Other Products We Offer

- **Sterilizer Charts:** Getinge & Steris
- **Lab Pharmacy:** Steris, Panasonic, Fujitsu, Linear Instruments, Pyxis, Seiko
- **Facilities Maintenance Charts:** Honeywell, Weksler, Bailey, United Electric, Instron
- **Variety of Pens and Markers**

0084845C



0099041C





SONY

Sony ultrasound paper provides exceptional image detail, excellent resolution and longevity of the printed image. All of these outstanding qualities enable Sony to offer the most trusted products to explore what is happening inside the body with great precision; clinicians and patients always notice the difference.



MOST POPULAR

UP110 Series BLACK & WHITE



The best selling **UPP-110 series** offers the following advanced features: high water resistance, head-matching performance, gray-scale reproduction, anti-electrostatic layer, minimal curling, high humidity and heat resistance, superior print quality and excellent tearing properties.

Product #	Cross Ref.	Description	Size	Per Pack	Per Case
2021977C	UPP110HA	Superior Density Printing Paper FOR USE IN: UP-880, UP-890MD	110mm x 18m	240 Prints	10 Rolls
2019086C	UPP110HD	High Density Printing Paper FOR USE IN: UP-890MD, UP0895MD, UP-D895, UP-D898MD, UP-X898MD	110mm x 20m	270 Prints	10 Rolls
0098226C	UPP110HG	High Gloss Media FOR USE IN: UP-895MD, UP-D895, UP-D898MD, UP-X898MD	110mm x 18m	240 Prints	10 Rolls
2019085C	UPP110S/5	Standard Media FOR USE IN: UP-890MD, UP-895MD, UP-D895, UP-897MD, UP-D897, UP-D898MD, UP-X898MD	110mm x 20m	270 Prints	10 Rolls

ULTRASOUND FILMS

Vermed is an authorized distributor for Sony and Mitsubishi diagnostic imaging media, providing exceptional image detail, excellent resolution, and longevity. We also offer Vermed Ultra 110HD and Ultra 110HG as competitive alternatives compatible with Sony printers.

SONY - Black and White Ultrasound Films

Product #	Cross Ref.	Description
2026054C	UPT210BL	Blue Thermal Transparency Film FOR USE IN: UP-980, and UP-991AD
0104937C	UPP725	8 x 10 Print Paper FOR USE IN: UP-D72XR
0106885C	UPT735BL	8 x 10 Blue Transparency Film FOR USE IN: UP-D71XR
7G31422960	UPP210SE	Standard Black & White Media FOR USE IN: UP-910, UP-930, UP-960, UP-980, UP-971AD, and UP-991AD
2026053C	UPP210HD	High Density Black & White Media FOR USE IN: UP-910, UP-930, UP-960, UP-980, UP-971AD, and UP-991AD

SONY - Colour Ultrasound Films

Product #	Cross Ref.	Description
7G32017149	UPC510C	A5 Colour Print Pack with Disposable Ribbon Holder FOR USE IN: UP-50, UP-D50, and UP-51MD
0104325C	UPC21L	A6 Colour Print Pack FOR USE IN: UP-21MD, UPD23MD, UP20, and UP25
0106244C	UPC21S	A6 Colour Print Pack FOR USE IN: UP-21MD, and UP-D21MD
7G32017030	UPC510	A5 Colour Print Pack FOR USE IN: UP-50, UP-D50, and P-51MD
7G32017150	UPC55	A5 Colour Print Pack FOR USE IN: UP-D55, and UP-D55MD

ULTRASOUND FILMS

MITSUBISHI

Mitsubishi ultrasound paper is a high quality product providing good resolution for ultrasound diagnostic printing. This paper provides great image detail when used with Mitsubishi printers such as: CP-700DSU, CP-700DU, CP-700UM, CP-770DW, P-68U, P-90W, P-91DW, P93DW, P-93W and P95DW.



7G30751336

Mitsubishi - Monochrome Ultrasound Films

Product #	Cross Ref.	Description	Size	Per Pack	Per Case
2018832C	K-61B	Super Grade Thermal Paper FOR USE IN: P-40U, P-50U, P-63UM, P-65U, P-67UA, P-68U, P-90U, P-91DW, P-91W, P93W, and P-95DW Mitsubishi Printers.	4.13" x 5.83"	100 Prints	20 Rolls
7G30957073	KP-65HM	High Density Thermal Sensitive Paper - Blue FOR USE IN: CP-700DSU, CP-700U, CP-700UM, CP-770DW, P-68U, P-90W, P-91DW, P-91W, P-93DW, P-93W, and P-95DW Mitsubishi printers.	110mm x 20m	215 Prints	20 Rolls
2020731C	K-75HM	High Density Thermal Paper - Blue FOR USE IN: P-71U, P-78U, and P-500W Mitsubishi printers	220mm x 25m	115 Prints	20 Rolls
0106115C	K-91HG	High Gloss Thermal Paper FOR USE IN: P-91DW, P93DW, P-93W and P-95DW Mitsubishi Printers.	110mm x 20m		20 Rolls

Mitsubishi - Colour Ultrasound Films

Product #	Cross Ref.	Description	Size	Per Pack	Per Case
2023136C	CK-700S	-	-	-	-
7G32017137	CK-700	Colour Roll Paper. 10mm x 22m FOR USE IN: CP-700DSU, CP-700DU, CP-700UM, CP-710U, and CP-770DW Colour printers	4" x 6"	-	10 Rolls
4008756C	CK-700L4P	Colour 4 Panel Paper Roll and Ink Ribbon FOR USE IN: CP-700DSU, CP-700DU, CP-700UM, CP-710U, CP-750U, and CP-770DW Colour printers	5.16 x 3.75"	75 Sheets	10 Rolls
0105985C	CK-3020L4P	-	-	-	-
2023137C	CK-700L	3 Panel Colour Paper Roll and Ink Sheet FOR USE IN: CP-800UM Colour printers	4 1/2 x 5"	160 Prints	10 Packs
4006030C	CK-800L	3 Panel Colour Paper Roll and Ink Sheet FOR USE IN: CP-800UM Mitsubishi printers	5" x 7"	100 Prints	10 Packs
0106262C	CK-900L	Standard L Size Ink Sheet FOR USE IN: CP-900DW, CP-900UM, and CP-910U Colour printers.	4.3" x 6.2"	130 Prints	10 Packs
0106263C	CK-900S4P	-	-	-	-
0104170C	CK-900S	Standard S Size Ink Sheet FOR USE IN: CP-900DW, CP900UM, and CP910U Colour printers.	4.1" x 4.3"	200 Sheets	10 Packs
4007230C	CK-3812	-	-	-	-
2019339C	CK-50S	-	-	-	-

ULTRASOUND FILMS

VERMED ULTRA

Available in both high density and high gloss formats, Vermed **Ultra 110HD** and **Ultra 110HG** are competitively priced thermal ultrasound films compatible with Sony medical grade printers.


Ultra 110HD - High Density Ultrasound Film

Product #	Rolls / Box	Rolls / Case
7G32019065	5 Rolls	50 Rolls

Ultra 100HG - High Gloss Ultrasound Film

Product #	Rolls / Box	Rolls / Case
7G32019066	5 Rolls	50 Rolls





Disposable CTG Belts: Disposable belts come in pairs of 2, 1 pink and 1 blue, with 2 buttons included.

Reusable CTG Belts: Reusable belts come in pairs of 2 Black Belts. Buttons are sold separately.

CTG / Abdominal Transducer Belts

Disposable and reusable fetal belts for secure CTG monitoring.

Product #	Cross Ref.	Description	Per Box
4009043C	DISPBELTBUT02	Disposable Abdominal Transducer Belt, 2 Pack - 1 x Pink, 1 x Blue, with 2 Buttons	100 Belts and 100 Buttons / 50 Packaged Pairs
4001531C	ROLLCTGBELT02	Disposable Abdominal Transducer Belt on a roll 52mm x 15m Black	-
4001530C	CTGBELT02	Reusable Abdominal Transducer Belt, 2 Pack Black	200 Belts / Box 100 Pairs

FETAL BELTS

Designed for ease of use for clinicians and comfort for patients, these Transducer Belts have easy lock/unlock button holes, placed at intervals to obtain the most comfortable and safely secured CTG monitor position.

4001531C

4001530C

4009043C



Vermed

Specially formulated to ensure excellent signal transmission.

Product #	Cross Ref.	Description
4009235C (Pictured Right)	G250B	250 ml Bottle - Blue ultrasound Gel



Aquasonic 100

The most widely used gel for diagnostic and therapeutic medical ultrasound procedures.

Product #	Cross Ref.	Description
4003889C (Pictured Right)	01-50	5L Sonicpac w/ Refillable Dispenser
4003891C	01-08	0.25L Dispenser
4003892C	01-02	60g Doppler sized tube
4007229C	01-34	1L w/ Dispenser and Cap



Aquasonic Clear

Contains no fragrance or colour. For diagnostic and therapeutic medical ultrasound procedures

Product #	Cross Ref.	Description
4000516C (Pictured Right)	03-50	5L Sonicpac w/ Refillable Dispenser
4008448C	03-08	0.25L Dispenser
4009439C	03-02	60g Doppler sized tube



ULTRASOUND GELS

Vermed carries an extensive line of ultrasound gels and accessories. From the most widely used gel for diagnostic ultrasound - Aquasonic 100, to Gel Warmers and Disinfectants. We carry all the gels and accessories required throughout your facility.



4000516C

4003889C

4009235C

SCAN Ultrasound Gels

Less viscous gel.

Product #	Cross Ref.	Description
4009440C	11-08	237ml (8 fl oz.) Dispenser
4009329C (Pictured Right)	11-28	3.8L (1 Gallon) w/ Refillable Dispenser



Ultrasound Probe Covers

Fits most endocavity ultrasound probes.

Product #	Cross Ref.	Description
4006729C	38-01	Latex Free Probe Cover
7G32021026	38-03	3D/4D Probe Cover



Thermasonic Gel Warmers

Single and 3 Bottle Warmers.

Product #	Cross Ref.	Description
7G32021021	82-01-20 CE	Single Bottle Warmer - 230V (Euro Plug)
7G32021022	82-03-20 CE	Multi Bottle Warmer - 230V (Euro Plug)
7G32024577 (Pictured Right)	82-03-20 CE	Multi Bottle Warmer LCD - 230V (Euro Plug)




Electromedical Gels, Pastes and Adhesives

Product #	Cross Ref.	Description
400920 (Pictured Right)	12-02	Spectra 360 Electrode Gel - 60g (2oz tube)
2000764C (Pictured Right)	12-08	Spectra 360 Electrode Gel - 250g (8.5oz tube)
7G32021031	22-06	Tensive Conductive Adhesive - 50g tube
7G32020230	67-05	Redux Electrolytic Paste - 5oz tube
2026134C	14-25	Signagel Electrode Gel - 250g (8.5oz tube)



A WEBSITE JUST FOR YOU

For more information and specifications visit us at Vermed.co.uk



COME CHECK US OUT

Visit our website to access all of our information. Search for products, download spec sheets and product images, and find digital copies of our catalogs. Our website is the place to go for any product information.

Visit us today at Vermed.co.uk to start your Vermed experience and feel confident that you're backed by the leading Sales & Customer Service Team.





Vermed, a Graphic Controls company is the ORIGINAL chart paper manufacturer and most diverse producer of electrodes for all departments in the hospital. Vermed manufactures in the USA and around the world to strict Original Equipment Manufacturer specifications and design criteria.

Welcome to the Vermed Experience

Vermed, a Graphic Controls company operates ISO 9001 and ISO 13485 Quality registered FDA registered facilities.

The OEM trademarks identified herein are the trademarks of the respective OEMs, and not Vermed. Vermed disclaims any affiliation, connection or association between it's products and those of the respective OEMs.



GLOBAL HEADQUARTERS
400 Exchange Street
Buffalo, NY 14204
United States
graphiccontrols.com
1.800.669.6905



BELLOWS FALLS, VERMONT
9 Lovell Drive
Bellows Falls, VT 05101
United States
vermed.com
1.800.245.4025



TOTNES, UNITED KINGDOM
St. Peter's Quay
Totnes, Devon TQ9 5XH
United Kingdom
vermed.co.uk
+44 1803 860100



INNSBRUCK, AUSTRIA
Egger-Lienz-Str. 1 D
6020 Innsbruck
Austria
vermed.co.uk
+43 512 94034



GENT, BELGIUM
Pijphoekstraat 4
B-9041 Gent(Oostakker)
Belgium
vermed.co.uk
+32 9228 1046



SANDPOINT, IDAHO
814 Airport Way
Sandpoint, Idaho 83864
United States
biomedical-innovations.com
1.800.201.3958



GANANOQUE, CANADA
23 Mill Street, Suite 201
Gananoque, Ontario K7G2L5
Canada
graphiccontrolsmedical.com
1.800.267.9498



DÜREN, GERMANY
Industriestrasse 8
D-52355 Düren (Lendersdorf)
Germany
vermed.co.uk
+49 / (0) 24 21 59 01 - 0



WARSAW, POLAND
ul. Sienkiewicza 22
05-420 Józefów
Poland
vermed.co.uk
+48 22 7895339

